

BrdU-ISH double labeling method

Akiya Watakabe (watakabe@nibb.ac.jp)

Here I wrote non-fluorescent and fluorescent protocols to double-label BrdU and ISH. These methods are to be combined with our protocol for normal ISH. BrdU is a Thymidine analogue and is now widely used to label dividing cell populations. For example, in the field of developmental neuroscience, it has been used in place of 3H-thymidine for neuronal birthdating. That is, particular types of neurons are often generated at certain developmental timepoints and stop dividing, and thus can be labeled by BrdU by simply injecting the pregnant animals with BrdU solution at those timepoints. BrdU labeling is also used to identify the newly generated neurons in the adult brains.

<Reagents>

BrdU solution (20 mg/ml) BrdU (sigma#B5002) in 0.007 N NaOH, 0.9% NaCl
PBST PBS with 0.2 % TritonX100
Anti-BrdU Abcam: BrdU antibody [BU1/75 (ICR1)]
Elite ABC kit Vectastain: PK-6100
Streptavidin-Cy3 Jackson ImmunoResearch laboratory: #016-160-084
TNT TS7.5, 0.05% Tween20

<BrdU labeling>

Inject BrdU solution into pregnant animals (50 mg/kg, i.p.).
Fix the offsprings at desired timepoints with 4% paraformaldehyde/0.1M PB by cardiac perfusion. Dissect out the brain and cut sections.

BrdU-ISH method (DAB-NBT/BCIP)

Fix in 4% PFA /0.1M PB at 4°C o/n

0.1M PB 10min x2

0.75% Glycine/0.1M PB 15min x2

0.3% Triton X100/0.1M PB 20min

0.1M PB 5min

ProK at 37°C 30min

acetylation 10min

0.1M PB 10min x2

Prehybridization at 60°C 1hr

Hybridization at 60°C o/n

Wash 2xSSC/FA/NLS at 60°C 15min x2

RNaseA Buffer 5min

RNaseA(20ug/ml) at 37°C 30min

2xSSC/NLS at 37°C 15min x2

0.2xSSC/NLS at 37°C 15min x2

TS7.5

Blocking 1%Blocking Reagent/TS7.5 30min

anti-DIG-AP (1/1000) in 1%Blocking at RT 2hr or at 4°C o/n

Wash TNT 15min x3

TS9.5 15min

NBT/BCIP(1/50) in TS9.5 1hr

Wash PBS 10min x2

1.5N HCl at 37°C 30min

Wash TNT 5min x2

Blocking 5% skim milk in PBST(SM/PBST) 20min

anti-BrdU(1/75) in SM/PBST at 4°C o/n

Wash PBS 10min x3

anti-rat-biotin(1/500) in SM/PBST at RT 2hr or at 4°C o/n

Wash TNT 10min x3

ABC(A:1/100, B:1/100, in TNT) 30min

Wash TNT 10min x3

DAB(H₂O₂:1/1000, DAB:1/100, in TNT) 10min

Wash PBS

BrdU-ISH method (Fluorescent)

Fix in 4% PFA /0.1M PB at 4°C o/n

0.1M PB 10min x2

0.75% Glycine/0.1M PB 15min x2

0.3% Triton X100/0.1M PB 20min

0.1M PB 5min

ProK at 37°C 30min

acetylation 10min

0.1M PB 10min x2

Prehybridization at 60°C 1hr

Hybridization at 60°C o/n

Wash 2xSSC/FA/NLS at 60°C 15min x2

RNaseA Buffer 5min

RNaseA(20ug/ml) at 37°C 30min

2xSSC/NLS at 37°C 15min x2

0.2xSSC/NLS at 37°C 15min x2

TS7.5

Blocking 1%Blocking Reagent/TS7.5 30min

anti-DIG-HRP (1/2000) in 1%Blocking at RT 2hr or at 4°C o/n

Wash TNT 15min x3

TSA-plus(1/50) 30min

Wash TNT 10min x3

1.5N HCl at 37°C 30min

Wash TNT 5min x2

Blocking 5% skim milk in PBST(SM/PBST) 20min

anti-BrdU(1/75) in SM/PBST at 4°C o/n

Wash PBS 10min x3

anti-rat-biotin(1/500) in SM/PBST at RT 2hr or at 4°C o/n

Wash TNT 10min x3

ABC(A:1/100, B:1/100, in TNT) 30min

Wash TNT 10min x3

anti-streptavidin-Cy3(1/500), anti-DNP-Alexa488(1/500) in SM/PBST at RT 2hr or at 4°C o/n

Wash TNT x3