

## Cryopreservation of *Chlamydomonas reinhardtii*

Ryutaro Tokutsu and Jun Minagawa

Ver. 0.9

Institute of Low Temperature Science, Hokkaido University  
North 19 West 8, Sapporo 060-0819, JAPAN  
Email: minagawa@lowtem.hokudai.ac.jp

### Freezing :

1. Cells are grown in a liquid TAP medium ( $1-6 \times 10^6$  cells / mL)
2. Mix 4.85 mL of culture and 150 mL of methanol  
(final conc. 3%) in a test tube
3. Transfer 500 mL of the mixture to a Falcon cryotube
4. Place the Falcon tube onto a "Mr. Frosty" (ambient temp.)
5. Place the "Mr. Frosty" in a deep freezer ( $-80^{\circ}\text{C}$ ) for 90 min
6. Cool down the Falcon tube to  $-196^{\circ}\text{C}$  in a liquid  $\text{N}_2$  bath
7. Store the frozen Falcon tubes in a liquid  $\text{N}_2$  container (Cryosystem 4000, MVE)

### Thawing :

8. Thaw a frozen stock at  $35^{\circ}\text{C}$
9. Transfer the thawed cell suspension into a TAP liquid medium (10mL) in a test tube
10. Incubate at  $25^{\circ}\text{C}$  over night under the dim light
11. Spread 800 mL of the culture on a TAP agar plate
12. Incubate at  $25^{\circ}\text{C}$  for 1-2 weeks
13. Typically 2000-4000 colonies appear within 1-2 weeks