

The Objectives of the Nikon A1R in Myodaiji area

No.	Name	Mag.	Immersion	N. A.	W.D. (mm)
1	PlanApo 10x	10	None	0.45	4.00
2	PlanApo-VC 20x	20	None	0.75	1.00
3	LWD Lambda S 40xWI	40	Water	1.15	0.61- 0.59
4	PlanApo-VC 60xA WI	60	Water	1.20	0.31 - 0.28
5	PlanApo-VC 60xH	60	Oil	1.40	0.13
6	PlanApo-VC 100xH	100	Oil	1.40	0.13

Mag.: Magnification, N. A.: Numerical Aperture, W. D.: Working Distance

The additional objectives for Nikon A1R

No.	Name	Mag.	Immersion	N. A.	W.D. (mm)
7	PlanApo 4x	4	None	0.2	20
8	PlanFluor 20xMI	20	Oil/Water/Glycerol	0.75	0.35
9	PlanApo 40x	40	None	0.95	0.16 - 0.12

Note: These objectives are kept in the lens cabinet in the microscope room. Please contact SBI, if you want to use. (Don't change it by yourselves.)