

The filter for visual observation of the Nikon A1R in Myodaiji area

Name	Excitation filter (nm)	Dichroic mirror (nm)	Emission filter (nm)
DAPI	340-380	400	435-485
FITC-LP *	465-495	505	510
TxRed	540-580	595	600-660
CFP	436/20	455	480/40
DAPI	340-380	400	435-485

The additional filter for the Nikon A1

Name	Excitation filter (nm)	Dichroic mirror (nm)	Emission filter (nm)
FITC-BP *	465-495	505	515-555

This filter cube is kept in the lens cabinet in the microscope room. Please contact SBI, if you want to use it.